

SET0125



CT compatible and MR conditional basic brachytherapy tool for use in various body sites. Well proven and easy to use.

Vagina

Ease of Use

The CT/MR Vaginal Cylinder Applicator Set contains the classical and proven vaginal applicators for intracavitary brachytherapy. Simple, yet effective, this applicator set can be used for gynecological treatments.

Easy Insertion

The cylinders are designed for easy insertion, and are tolerable for the patient during treatment. Cylinders are available in four different diameters to easily adapt to the individual patient's anatomy. The closed and curved end of the applicator facilitates insertion into the body.

Precise Dose Distribution

The CT compatible and MR conditional applicator set allows state of the art planning, improving the treatment planning processes. The cylinder design provides the desired distances between the source and target area. This applicator set is an essential brachytherapy tool for gynecological treatments.

Features

- CT compatible and MR conditional
- Cylinders of different diameters adapt to patient anatomy
- Design facilitates insertion
- Defined geometry of cylinders ensures homogeneous dose distribution at the surface
- Set is reusable; reprocessing via automatic cleaning and disinfection; steam sterilization
- Applicator can easily be attached to the applicator clamp (LAZ22-01)

Set Content SET0125

Item	Part number
Vaginal Cylinder, closed, D=22 mm, L=130 mm	LA320-64
Vaginal Cylinder, closed, D=25 mm, L=130 mm	LA320-74
Vaginal Cylinder, closed, D=30 mm, L=130 mm	LA320-84
Vaginal Cylinder, closed, D=35 mm, L=130 mm	LA320-89
Adapter for Vaginal Cylinder	LC313-22
Fixation Tube, D=5 mm	LC313-21


Required Items

Alternative 1	Part number
Transfer Tube, D = 3 mm, L = 1000 mm	LAF1000
CT/MR Intrauterine Tube, straight, thin, D=3.5 mm	LCT 01-01
Alternative 2	Part number
Universal Applicator incl. plug, D=3 mm, L=1400 mm (non sterile)	LAA1400-GYN
Intrauterine Tube, CT/MR, straight	LCR01-01

Optional Items

Item	Part number
Sterilization Box, 381 x 254 x 38 mm ¹	LLZ05-01

Isodose



Transfer Tube



Sterilization Box



¹ Manufactured by PST Plastic Sterilizing Trays Corporation

* For complete information and prior to use please read the Instructions for Use

Manufacturer:

BEBIG Medical GmbH

Robert-Rössle-Str. 10
13125 Berlin
Germany

Phone +49 30 9487888-10

sales@bebigmedical.com

Regional Sales, Marketing and Services:

Global

BEBIG Medical GmbH

Robert-Rössle-Str. 10
13125 Berlin
Germany

Phone +49 30 9487888-10

sales@bebigmedical.com

North America

Mick Radio-Nuclear Instruments, Inc. A BEBIG Medical Company

521 Homestead Avenue
Mount Vernon, NY 10550
USA

Phone +1 914 667 3999

sales@micknuclear.com

CE 0344

www.bebigmedical.com